Engenharia

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
(J. R. Drummond, 1996)

General Editing

Pop-Up	Right-click
Rotate Tools	tab
Operating ⇌ Positioning (front panel) Positioning ⇌ Wiring (block diagram)	space
Move item one (several) pixel(s)	(Shift)-arrow
Increase/Decrease Controls faster	shift-click
add to ring control	shift-enter
Toggle help on/off	Ctrl-h
duplicate using positioning Tool	Ctrl-drag

Diagram Editing

Create subVI when repeating part of a block diagram	Select items Edit>> Create SubVI from Selection
Create a constant/control/indicator	Pop-up on function terminals Create Constant/Control/Indicator
Create a constant	drag a front panel control to block diagram
Create a control	drag block diagram constant to front panel
List errors	Pop-up on bad wires, List Errors -or- click on broken arrow icon
Open VI front panel (block diagram)	(Ctrl)-Double-click
Delete bad wires	Ctrl-b


Help

Show Help	Context-Sensitive (Ctrl-h)
Lock Help	Hold that help! (Ctrl-Shift0-L)
Simple Help	Controls amount of help shown
Online Reference	Complete documentation (Ctrl-?)
Online Help	Customised Online reference for this VI
About LabVIEW	Software rev. & serial #


Wiring Techniques

Tack down a wire	Click
Un-tack last point	Ctrl-Click
Switch wire direction	Space
Highlight segment/wire/net	Click/Double-click/triple-click
Tack end of wire (temporary)	Double-click


Tools

Controls

Functions

LabVIEW Documentation and Help Resources

Context Help

When you move the cursor over objects on the front panel or block diagram, the **Context Help** window displays basic information about each object. Select **Help»Show Context Help** to display the **Context Help** window.

LabVIEW Help

Includes reference information for each palette, menu, tool, dialog box, and built-in VI or function. The *LabVIEW Help* also includes step-by-step instructions for using LabVIEW features. Select **Help»VI, Function, & How-To Help** to access the *LabVIEW Help*.

LabVIEW Bookshelf

Contains links to all the LabVIEW documentation in electronic format. Use the *LabVIEW Bookshelf* to search PDF versions of all the LabVIEW manuals and Application Notes. Select **Help»Search the LabVIEW Bookshelf** to access the *LabVIEW Bookshelf*.

LabVIEW Examples

LabVIEW includes hundreds of example VIs you can use and incorporate into your own VIs. You can modify an example to fit your application, or you can copy and paste from one or more examples into your own VI. Browse or search the example VIs by selecting **Help»Find Examples**.

LabVIEW Web Resources

Includes a list of links to online National Instruments resources. Select **Help»Web Resources** to access the LabVIEW resources at ni.com. Resources available to you include the following:

Products and Services

- Upgrades
- Training courses
- Seminars and user groups

Solutions

- Systems integrators
- Customer solutions

Support

- LabVIEW support
- Manuals

NI Developer Zone

- Discussion forums
- Example code

ni.com/labviewzone

Keyboard Shortcuts

Refer to the LabVIEW Help to access this list of keyboard shortcuts.

Object/Movement

Shift-click	Selects multiple objects; adds object to current selection.
N++ (arrow keys)	Moves selected objects one pixel at a time.
Shift-N++	Moves selected objects several pixels at a time.
Shift-click (drag)	Moves selected objects in one axis.
Ctrl-click (drag)	Duplicates selected object.
Ctrl-Shift-click (drag)	Duplicates selected object and moves it in one axis.
Shift-resize	Resizes selected object while maintaining aspect ratio.
Ctrl-resize	Resizes selected object while maintaining center point.
Ctrl-drag a rectangle	Inserts more working space on the front panel or block diagram.
Ctrl-A	Performs last alignment operation on objects.
Ctrl-D	Performs last distribution operation on objects.
Double-click open space	Places free label on the front panel or block diagram if automatic tool selection is enabled.
Ctrl-mouse wheel	Scrolls through subdiagrams of a Case, Event, and Stacked Sequence structure.

Navigating the LabVIEW Environment

Ctrl-E	Displays block diagram or front panel.
Ctrl-G	Enables or disables alignment grid. On French keyboards, press the <Ctrl-;> keys.
Ctrl-J	Maximizes and restores window.
Ctrl-T	Tiles front panel and block diagram windows.
Ctrl-F	Finds objects or text.
Ctrl-G	Searches VIs for next instance of object or text.
Ctrl-Shift-G	Searches VIs for previous instance of object or text.
Ctrl-Shift-F	Displays Search Results window.
Ctrl-Tab	Cycles through LabVIEW windows.
Ctrl-Shift-Tab	Cycles backward through LabVIEW windows.
Ctrl-I	Displays VI Properties dialog box.
Ctrl-L	Displays Error Use window.
Ctrl-Y	Displays History window.
Ctrl-D	Redraws Hierarchy window.
Ctrl-A	Shows all VIs in Hierarchy window.
Enter *	Finds next node that matches the search string.
Shift-Enter *	Finds previous node that matches the search string.

* After initiating a search by typing in Hierarchy window

Debugging

Ctrl-I	Steps into node.
Ctrl-→	Steps over node.
Ctrl-t	Steps out of node.

File Operations

Ctrl-N	Creates new VI.
Ctrl-O	Opens existing VI.
Ctrl-W	Closes VI.
Ctrl-S	Saves VI.
Ctrl-P	Prints the window.
Ctrl-Q	Quits LabVIEW.

Basic Editing

Ctrl-Z	Undoes last action.
Ctrl-Shift-Z	Redoes last action.
Ctrl-X	Cuts an object.
Ctrl-C	Copies an object.
Ctrl-V	Pastes an object.

Help

Ctrl-H	Displays Context Help window.
Ctrl-Shift-L	Locks Context Help window.
Ctrl-? or F1	Displays LabVIEW Help.

Note The <Ctrl> key in these shortcuts corresponds to the (Mac OS) <Option> or <Command> key, (Sun) <Meta> key, or (Linux) <Alt> key.

Tools and Palettes

Ctrl	Switches to next most useful tool.
Shift	Switches to Positioning tool.
Ctrl-Shift over open space	Switches to Scrolling tool.
Spacebar *	Toggles between two most common tools.
Shift-Tab *	Enables automatic tool selection.
Tab *	Cycles through tools if you disabled automatic tool selection by clicking the Automatic Tool Selection button. Otherwise, enables automatic tool selection.
F4++	Navigates temporary Controls and Functions palettes.
Enter	Navigates into temporary palette.
Esc	Navigates out of temporary palette.

* If automatic tool selection is disabled

SubVIs

Double-click subVI	Displays subVI front panel.
Ctrl-double-click subVI	Displays subVI block diagram.
Drag VI icon to block diagram	Places that VI as a subVI on the block diagram.
Shift-drag VI icon to block diagram	Places that VI as a subVI on the block diagram with constants wired for controls that have non-default values.
Ctrl-right-click block diagram and select VI	Opens front panel of that VI from the Functions palette.

Execution

Ctrl-R	Runs VI.
Ctrl-. *	Aborts VI.
Ctrl-M	Changes to run or edit mode.
Ctrl-Run button	Recompiles current VI.
Ctrl-Shift-Run button	Recompiles all VIs in memory.
Ctrl-I *	Moves key focus inside an array or cluster.
Ctrl-O *	Moves key focus outside an array or cluster.
Tab *	Navigates the controls or indicators according to tabbing order.
Shift-Tab *	Navigates backward through the controls or indicators.

* While VI is running

Wiring


Ctrl-B	Removes all broken wires.
Esc, right-click, or click terminal	Cancels a wire you started.
Single-click wire	Selects one segment.
Double-click wire	Selects a branch.
Triple-click wire	Selects entire wire.
A	Disables automatic wire routing temporarily.
Double-click (while wiring) spacebar	Tacks down wire without connecting it.
Shift-click	Toggles automatic wiring while moving objects.
Ctrl-click	Undoes last point where you set a wire.
Ctrl-click Input on function with two inputs	Switches the two Input wires.
spacebar	Switches the direction of a wire between horizontal and vertical.

Text

Double-click	Selects single word in string.
Triple-click	Selects entire string.
Ctrl-→	Moves forward in string by one word.
Ctrl-←	Moves backward in string by one word.
Home	Moves to beginning of current line in string.
End	Moves to end of current line in string.
Ctrl-Home	Moves to beginning of entire string.
Ctrl-End	Moves to end of entire string.
Shift-Enter	Adds new items when entering items in enumerated type controls and constants, ring controls and constants, or Case structures.
Esc	Cancels current edit in a string.
Ctrl-Enter	Ends text entry.
Ctrl=	Increases current font size.
Ctrl-	Decreases current font size.
Ctrl-F	Displays Font dialog box.
Ctrl-1 *	Displays Application font.
Ctrl-2 *	Displays System font.
Ctrl-3 *	Displays Dialog font.
Ctrl-4 *	Displays current font.

* In Font dialog box

Editing, Execution, and Debugging Tools

 Show Context Help Window —Displays the Context Help window.	 Resize Objects —Resizes multiple front panel objects to the same size.
 Run —Runs the VI.	 Reorder —Reorders objects if they overlap, including moving forward and moving backward.
 Broken Run —Indicates that the VI contains errors. Click the button to list errors.	 Highlight Execution —Shows the flow of data through the block diagram.
 Run Continuously —Runs the VI until you abort or pause it.	 Step Into —Steps into a node, which blinks when ready to execute.
 Abort Execution —Stops the VI immediately.	 Step Over —Executes a node and pauses at the next node.
 Pause —Pauses or resumes a running VI.	 Step Out —Steps out of a node and continues to the next node.
 Align Objects —Aligns objects along axes.	 Breakpoint Tool —Sets breakpoints on nodes to pause execution.
 Distribute Objects —Spaces objects evenly.	 Probe Tool —Creates probes on wires to display intermediate values in a running VI.
 Text Settings —Changes the attributes of text in the VI.	 Call Chain —Lists all the callers from the top-level VI to the opened subVI.

Terminal Data Types

Signed Integers (i)	8-bit 	Boolean (FALSE) 	 Waveform
	16-bit 	String (empty string) 	 Digital waveform
	32-bit 	Path (<Not A Path>) 	 Digital
Unsigned Integers (u)	8-bit 	Reference number 	 I/O name
	16-bit 	Enumerated type 	 Time stamp
	32-bit 	Cluster	Numeric 
Floating-Point (f)	Single 	Mixed data 	 Dynamic
	Double 	Array	 P closure
	Extended 	1D 	 Variant
		2D 	
Complex Floating-Point (c + 0.0)	Single 		
	Double 		
	Extended 		

 Note: Default values appear in ().